

LUNATIC



LUNCHAINED
LABS

Get your quant on

Lunatic makes batch quantification of protein, DNA and RNA a no-brainer. All you need is 2 μ L and 10 minutes to measure up to 96 samples. Run them straight-up, even at high concentrations, without ever diluting. Lunatic gets biologics and genomics UV/Vis quantification on the money every time. Just drop, load and read.



2 μ L volume


96 samples

10 minutes

SBS plate format



No B.S. workflow

Skip the sample prep, cleaning and worrying about cross-contamination or evaporation – your samples will sit tight for up to 2 hours. Our Lunatic plate runs 96 samples in 10 minutes. Each microfluidic circuit has two fixed pathlengths built-in to cover a wide dynamic range of 0.03–275 OD. All you do is load up your samples and get kick-ass results.




Max out your biologic

Lunatic is the only system out there that can measure biologics at high-throughput and high concentration. It's got dynamic range that covers from 0.02 mg/mL to 200 mg/mL (mAb), so run any protein without ever having to dilute again. Stop the madness of running just one protein at a time – and cleaning up afterwards.



Ridiculously good data

Get spot-on precision within 1% and accuracy within 2%. Using two fixed pathlengths, Lunatic gets you jaw-dropping data at both low and high protein concentrations and nails the expected concentrations every time. With the smallest sample size, the highest throughput and crazy accurate results, Lunatic is hands down the best tool out there for protein quantification.




Known Conc. (mg/mL)	Average Conc. (mg/mL)	CV (%)
201.4	201.7	0.7%
169.9	171.1	0.7%
139.4	140.7	0.4%
113.1	115.0	0.6%
79.4	80.7	0.2%
50.2	51.2	0.2%
19.9	20.3	0.3%
12.1	12.2	0.2%
7.21	7.28	0.2%
2.35	2.37	0.2%

Accuracy you can show off

Lunatic has shocking accuracy and wants to prove it. Using the Fundamentals tryptophan standards, the accuracy, precision and linearity of the instrument can be proven at any time – at the protein-relevant 280 nm wavelength and OD range of 20-225. Be confident about every sample with data that leaves no room for doubt.




The Fundamentals
Certified tryptophan standards
OD range 20-225




Downstream ready

Smash through the USP & Ph. Eur. UV/Vis requirements for absorbance accuracy, precision, linearity, wavelength accuracy, stray light and resolution with performance verification measurements of independently certified NIST standards. Lunatic is ready when you are to get validated and make the move into QC.



Performance verification plate

Holds certified pharmacopeia standards for UV/Vis verification



Put it on lockdown

GLP labs don't sweat it. Lunatic's software hooks labs up with 21 CFR Part 11 compliant features. We're talking password protection, electronic signatures, full audit trail – the whole package.

The screenshot displays the 'Lunatic & Stunner Client' software interface. The top navigation bar is green and contains the application name, the identifier '21CFRp11', and the user role 'Administrator' with a profile icon. Below the navigation bar, there are two main sections: 'Settings' and 'Audit Trail'.

Settings

Force change password after (days): 60
Maximum inactivity time (minutes): 30
Maximum number of failed logins: 5
Minimum characters in password: 5

Password needs to contain:

- numbers
- upper case letters
- lower case letters
- special characters

A 'SAVE CHANGES' button with a floppy disk icon is located at the bottom right of the settings panel.

Audit Trail


Category: All | Action: All | EXPERIMENT'S INTEGRITY | SAVE TO PDF

ID	Date	User	Category	Action	Details
2019	17/04/2020 10:45:49	Administrator	User	User login	User login successful (Administrator)
2018	17/04/2020 10:45:34		System	Software start	Software start: Lunatic & Stunner Client
2017	17/04/2020 10:45:10		System	Software shutdown	Software shutdown: Lunatic & Stunner Client
2016	17/04/2020 10:44:57	Administrator	System	Security settings modified	"Force change password after" changed from "90 days" to "60 days"
2015	17/04/2020 10:44:28	Administrator	System	CRM info modified	Absorbance CRM certified absorbances modified
2014	17/04/2020 10:43:59	Administrator	System	PV criteria modified	"Absorbance Linearity minimum R square" changed from "0.000" to "0.9
2013	17/04/2020 10:43:31	Administrator	System	CRM info modified	Absorbance CRM info concentration 1 modified: "serial number" change
2012	17/04/2020 10:41:03	Administrator	System	Add license	Seat license added (hardwareid: "*****A01214") which activates 21C
2011	17/04/2020 10:38:14	Administrator	User	User login	User login successful (Administrator)
2010	17/04/2020 10:37:16		System	Software start	Software start: Lunatic & Stunner Client
2009	17/04/2020 00:35:09		System	Software shutdown	Software shutdown: Lunatic & Stunner Analysis
2008	16/04/2020 23:24:42	Administrator	Data	New experiment	New experiment done on Fri, 20 Mar 2020 11:37:53, performed on S/N 4
2007	16/04/2020 22:41:18	Administrator	User	User login	User login successful (Administrator)

A green 'BACK' button with a left-pointing arrow is located at the bottom left of the interface.

Dig up your genomic's dirt

DNA and RNA samples can be freaking messy. Lunatic analysis software lets you see annoying impurities that classic A260 measurements alone mistake as your DNA or RNA. Ditch the time-sucking, dye-based prep and know right away if your samples are good to go.



Set it and forget it

Lunatic talks to robots. Hook it up with any liquid handler system and run your measurements completely hands-off with our API setup and keep more of your precious sample for the rest of your workflow.

Automated sample loading



Going smaller

Little Lunatic is there for labs that have fewer samples and lower requirements. It uses a chip-based version of the Lunatic plate to get quick concentration measurements of up to 16 samples in 5 minutes. Researchers that just need a quick check of concentration can get what they need from this little guy.



2 μ L volume

16 samples

5 minutes

Touch screen

Specifications

Lunatic instrument specifications		
Dimensions	37 cm W x 46 cm D x 33 cm H; 21 kg	
Operating voltage	24 VDC, 30 W (max)	
Connection	USB, TCP/IP (Service, SiLA)	
Light source	Xenon flash lamp	
Detectors	UV/Vis polychromatic spectrophotometer	
Approval	CE, FCC, CSA	
Wavelength range	230-750 nm	
Wavelength accuracy	≤400 nm: ±1 nm; ≥400 nm: ±2 nm	
Spectral resolution	Better than 2 nm (toluene in hexane)	
Absorbance precision (1 cm quartz cuvette)	<1 OD: ±0.005 OD st dev	1-2 OD: ±0.5% CV
Absorbance accuracy (1 cm quartz cuvette)	<1 OD: ±0.01 OD	1-2 OD: ±1%
Lunatic plate specifications	High Lunatic Plate	Lunatic Plate
Application area	High concentration samples, typically proteins	Low concentration samples, typically nucleic acids
Sample retention time	Up to 2 hours	Up to 2 hours
Recommended sample volume	2 µL	2 µL
Pathlength(s)	0.1 mm & 0.7 mm path	0.5 mm path
Measurement time for full plate	~10 minutes	~5 minutes
Measurement range: OD 10 mm ng/µL dsDNA mg/mL ave protein mix	0.03-275 OD 10 mm 1.5-13750 ng/µL 0.03-275 mg/mL	0.03-40 OD 10 mm 1.5-2000 ng/µL 0.03-40 mg/mL
Absorbance precision (10 mm pathlength)	<1 OD: ±0.01 OD st dev 1-200 OD: ±1% CV	<1 OD: ±0.01 OD st dev >1 OD: ±1.5% CV
Absorbance accuracy (10 mm pathlength)	<1 OD: ±0.02 OD 1-200 OD: ±2%	<1 OD: ±0.02 OD >1 OD: ±4%
Dimensions	8.55 cm W x 12.8 cm D x 1.21 cm H	
Samples per plate	96 (12 x 8 microplate format)	



Unchained Labs

6870 Koll Center Parkway

Pleasanton, CA 94566

Phone: 1.925.587.9800

Toll-free: 1.800.815.6384

Email: info@unchainedlabs.com

© 2020 Unchained Labs. All rights reserved. The Unchained Labs logo, Lunatic and the Lunatic logo are trademarks and/or registered trademarks of Unchained Labs. All other brands or product names mentioned are trademarks owned by their respective organizations.