



**VELOX**<sup>™</sup>  
Solvent Waste System  
Anti Static material




# Velox Solvent waste systems

## Introduction:

Fluids in HPLC waste containers are overwhelmingly dangerous containing potent volatile organic compounds (VOC) posing potent risk to health of Lab workers & lab safety.

Zodiac came with an innovative solution of introducing Anti Static material for **VELOX Solvent waste system** with activated carbon filters not only creating an organised containment for HPLC solvent waste but also providing an efficient way for collection and disposal of volatile vapours.



## Advantages

- Filters system against harmful vapours
- Capillary fixations with PFA Fittings
- Several systems can be connected to 1 container
- Different size containers compatible to any space condition





90 Adaptor included for maximum storage flexibility



Normal cap is additionally provided for convenient transport of the wastes



Blind fittings are provided to block the remaining ports when not in use

*Zodiac Velox Solvent Waste System allows 99 percent of the volatile vapours to be trapped into its Powerful exhaust filters keeping the laboratory environment safe and clean , This system is more organised than traditional techniques protects against spills and leaks. This is a necessary addition to every lab for improved solvent waste handling and disposal*

### Item Description

### Part Number

Velox Solvent Waste System - 5Ltrs

ZLS.332.2419

Velox Solvent Waste System - 10Ltrs

ZLS.332-3413

Velox Solvent Waste System - 20Ltrs

ZLS.332-4413

